

A1082

Leader in Biomolecular Solutions for Life Science



## FAM3B Rabbit pAb

Catalog No.: A1082

1 Publications

### Basic Information

#### Observed MW

#### Calculated MW

26kDa

#### Category

Polyclonal Antibody

#### Applications

IHC-P,IF/ICC,ELISA

#### Cross-Reactivity

Human,Mouse,Rat

### Background

Involved in insulin secretion. Located in extracellular exosome.

### Recommended Dilutions

IHC-P	1:100 - 1:200
IF/ICC	1:100 - 1:200

### Immunogen Information

#### Gene ID

54097

#### Swiss Prot

P58499

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 30-235 of human FAM3B (NP\_478066.3).

#### Synonyms

2-21; ORF9; PANDER; PRED44; C21orf11; C21orf76; FAM3B

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification


#### Storage

Store at -20°C. Avoid freeze / thaw cycles.


Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data


---




Immunohistochemistry analysis of paraffin-embedded human colon using FAM3B antibody (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded human colon carcinoma using FAM3B antibody (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of paraffin-embedded human stomach using FAM3B antibody (A1082) at dilution of 1:100 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunofluorescence analysis of mouse pancreas cells using FAM3B antibody (A1082) at dilution of 1:100. Blue: DAPI for nuclear staining.