

A11869

Leader in Biomolecular Solutions for Life Science



ICAM-1/CD54 Mouse mAb

Catalog No.: A11869

1 Publications

Basic Information

Observed MW

95kDa

Calculated MW

58kDa

Category

SMab Recombinant Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a cell surface glycoprotein which is typically expressed on endothelial cells and cells of the immune system. It binds to integrins of type CD11a / CD18, or CD11b / CD18 and is also exploited by Rhinovirus as a receptor.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3383

Swiss Prot

P05362

Immunogen

Recombinant protein of human ICAM1

Synonyms

BB2; CD54; P3.58; ICAM-1/CD54

Contact



www.abclonal.com

Product Information

Source

Mouse

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of lysates from Raji cells, using ICAM-1/CD54 Mouse mAb (A11869).
Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.



Western blot analysis of various lysates using ICAM-1/CD54 Mouse mAb (A11869).
Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) (AS003) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.