

A12267

Leader in Biomolecular Solutions for Life Science



## KIAA1524 Rabbit pAb

Catalog No.: A12267

4 Publications

### Basic Information

#### Observed MW

102kDa

#### Calculated MW

102kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human, Mouse, Rat

### Background

Enables protein homodimerization activity. Predicted to act upstream of or within positive regulation of neural precursor cell proliferation and spermatogenesis. Located in cytosol and plasma membrane.

### Recommended Dilutions

WB 1:500 - 1:1000

### Immunogen Information

#### Gene ID

57650

#### Swiss Prot

Q8TCG1

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 681-905 of human KIAA1524 (NP\_065941.2).

#### Synonyms

p90; NOCIVA; KIAA1524

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification


#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

## Validation Data

---



Western blot analysis of lysates from NIH/3T3 cells, using KIAA1524 Rabbit pAb (A12267) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 20s.