

A12377

Leader in Biomolecular Solutions for Life Science



## TEKT5 Rabbit pAb

Catalog No.: A12377

### Basic Information

**Observed MW**

56kDa

**Calculated MW**

56kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Mouse, Rat

### Background

Predicted to be involved in cilium assembly and cilium movement involved in cell motility. Located in nucleus.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

146279

**Swiss Prot**

Q96M29

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TEKT5 (NP\_653275.1).

**Synonyms**

CT149; TEKT5

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

## Validation Data

---



Western blot analysis of lysates from rat testis, using TEK5 Rabbit pAb (A12377) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 10s.