

A12378

Leader in Biomolecular Solutions for Life Science



TEKT5 Rabbit pAb

Catalog No.: A12378

Basic Information

Observed MW

56kDa

Calculated MW

56kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to be involved in cilium assembly and cilium movement involved in cell motility.
Located in nucleus.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

146279

Swiss Prot

Q96M29

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human TEKT5 (NP_653275.1).

Synonyms

CT149; TEKT5

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using TEKT5 Rabbit pAb (A12378) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.