

A12470

Leader in Biomolecular Solutions for Life Science



NDP Rabbit pAb

Catalog No.: A12470

Basic Information

Observed MW

Calculated MW
15kDa

Category
Polyclonal Antibody

Applications
ELISA

Cross-Reactivity
Rat

Background

This gene encodes a secreted protein with a cystein-knot motif that activates the Wnt/beta-catenin pathway. The protein forms disulfide-linked oligomers in the extracellular matrix. Mutations in this gene result in Norrie disease and X-linked exudative vitreoretinopathy.

Recommended Dilutions

Immunogen Information

Gene ID	Swiss Prot
4693	Q00604

Immunogen
Recombinant fusion protein containing a sequence corresponding to amino acids 25-133 of human NDP (NP_000257.1).

Synonyms
ND; EVR2; FEVR; NDP

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.