

A12495

Leader in Biomolecular Solutions for Life Science



RBP4 Rabbit pAb

Catalog No.: A12495

Basic Information

Observed MW

19kDa

Calculated MW

21kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

This protein belongs to the lipocalin family and is the specific carrier for retinol (vitamin A alcohol) in the blood. It delivers retinol from the liver stores to the peripheral tissues. In plasma, the RBP-retinol complex interacts with transthyretin which prevents its loss by filtration through the kidney glomeruli. A deficiency of vitamin A blocks secretion of the binding protein posttranslationally and results in defective delivery and supply to the epidermal cells.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5950

Swiss Prot

P02753

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 19-201 of human RBP4 (NP_006735.2).

Synonyms

RDCCAS; MCOPCB10; RBP4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using RBP4 Antibody (A12495) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.