

A12621

Leader in Biomolecular Solutions for Life Science



Occludin Rabbit pAb

Catalog No.: A12621

15 Publications

Basic Information

Observed MW

62kDa

Calculated MW

59kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

This gene encodes an integral membrane protein that is required for cytokine-induced regulation of the tight junction paracellular permeability barrier. Mutations in this gene are thought to be a cause of band-like calcification with simplified gyration and polymicrogyria (BLC-PMG), an autosomal recessive neurologic disorder that is also known as pseudo-TORCH syndrome. Alternative splicing results in multiple transcript variants. A related pseudogene is present 1.5 Mb downstream on the q arm of chromosome 5.

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

100506658

Swiss Prot

Q16625

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 350-450 of human Occludin (NP_002529.1).

Synonyms

BLC-PMG; PTORCH1; PPP1R115; Occludin

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.05% proclin300, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of 293T, using Occludin antibody (A12621) at 1:2000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 5s.