

A12634

Leader in Biomolecular Solutions for Life Science



## EVC Rabbit pAb

Catalog No.: A12634

### Basic Information

#### Observed MW

Refer to Figures

#### Calculated MW

112kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Human

### Background

This gene encodes a protein containing a leucine zipper and a transmembrane domain. This gene has been implicated in both Ellis-van Creveld syndrome (EvC) and Weyers acrodermal dysostosis.

### Recommended Dilutions

WB 1:1000 - 1:2000

### Immunogen Information

#### Gene ID

2121

#### Swiss Prot

P57679

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 743-992 of human EVC (NP\_714928.1).

#### Synonyms

EVC1; EVCL; DWF-1; EVC

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.