

A12668

Leader in Biomolecular Solutions for Life Science



## ZC3H12A Rabbit pAb

Catalog No.: A12668

### Basic Information

#### Observed MW

Refer to Figures

#### Calculated MW

66kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Mouse

### Background

ZC3H12A is an MCP1 (CCL2; MIM 158105)-induced protein that acts as a transcriptional activator and causes cell death of cardiomyocytes, possibly via induction of genes associated with apoptosis.

### Recommended Dilutions

WB 1:1000 - 1:2000

### Immunogen Information

#### Gene ID

80149

#### Swiss Prot

Q5D1E8

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZC3H12A (NP\_079355.2).

#### Synonyms

Reg1; MCP1P; MCP1P1; MCP1P-1; dJ423B22.1; ZC3H12A

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.