

A12741

Leader in Biomolecular Solutions for Life Science



Os06g0643000 Rabbit pAb

Catalog No.: A12741

Basic Information

Observed MW

Refer to Figures

Calculated MW

79kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Oryza sativa

Background

GPA5 is a membrane-associated protein capable of forming homodimers and that it is specifically localized to DVs in developing endosperm.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4341636

Swiss Prot

Immunogen

Recombinant protein of Oryza sativa Os06g0643000

Synonyms

Os06g0643000

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.