

A12746

Leader in Biomolecular Solutions for Life Science



IFNA8 Rabbit pAb

Catalog No.: A12746

Basic Information

Observed MW

Refer to Figures

Calculated MW

22kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse, Rat

Background

Predicted to enable cytokine activity and type I interferon receptor binding activity. Predicted to be involved in several processes, including B cell activation; lymphocyte activation involved in immune response; and positive regulation of peptidyl-serine phosphorylation of STAT protein. Predicted to act upstream of or within defense response to virus. Predicted to be located in extracellular region. Predicted to be active in extracellular space.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3445

Swiss Prot

P32881

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 24-189 of human IFNA8 (NP_002161.2).

Synonyms

IFN-alphaB; IFNA8

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.