

A12765

Leader in Biomolecular Solutions for Life Science



FAM46A Rabbit pAb

Catalog No.: A12765

1 Publications

Basic Information

Observed MW

50kDa

Calculated MW

50kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Enables RNA binding activity. Predicted to be involved in mRNA stabilization. Predicted to act upstream of or within response to bacterium. Implicated in lung non-small cell carcinoma; osteoarthritis; and osteogenesis imperfecta type 18.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

55603

Swiss Prot

Q96IP4

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 180-400 of human FAM46A (NP_060103.2).

Synonyms

OI18; XTP11; FAM46A; C6orf37

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using FAM46A Rabbit pAb (A12765) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit (RM00021). Exposure time: 90s.