

# FEM1B Rabbit pAb

Catalog No.: A12803

## Basic Information

### Observed MW

Refer to Figures

### Calculated MW

70kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Mouse

## Background

This gene encodes an ankyrin repeat protein that belongs to the death receptor-associated family of proteins and plays a role in mediating apoptosis. The encoded protein is also thought to function in the replication stress-induced checkpoint signaling pathway via interaction with checkpoint kinase 1.

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

10116

### Swiss Prot

Q9UK73

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 368-627 of human FEM1B (NP\_056137.1).

### Synonyms

F1AA; F1A-ALPHA; FEM1-beta; FEM1B

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.