

A12854

Leader in Biomolecular Solutions for Life Science



THAP3 Rabbit pAb

Catalog No.: A12854

Basic Information

Observed MW

27kDa

Calculated MW

27kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable DNA binding activity and metal ion binding activity. Involved in positive regulation of transcription by RNA polymerase II.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

90326

Swiss Prot

Q8WTV1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-239 of human THAP3 (NP_001182682.1).

Synonyms

THAP3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using THAP3 Rabbit pAb (A12854) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.