

A12887

Leader in Biomolecular Solutions for Life Science



## CT45A3 Rabbit pAb

Catalog No.: A12887

### Basic Information

**Observed MW**

21kDa

**Calculated MW**

21kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Mouse, Rat

### Background

Predicted to be involved in snRNA 3'-end processing. Predicted to be part of integrator complex.

### Recommended Dilutions

**WB** 1:1000 - 1:2000

### Immunogen Information

**Gene ID**

441519

**Swiss Prot**

Q8NHU0

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-189 of human CT45A3 (NP\_001017435.1).

**Synonyms**

CT45-3; CT45-4; CT45.3; CT45.4; CT45A4; CT45A3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

## Validation Data

---



Western blot analysis of various lysates using CT45A3 Rabbit pAb (A12887) at 1:3000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 90s.