

A12907

Leader in Biomolecular Solutions for Life Science



BMF Rabbit pAb

Catalog No.: A12907

Basic Information

Observed MW

Calculated MW
21kDa

Category
Polyclonal Antibody

Applications
ELISA

Cross-Reactivity
Human

Background

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. This protein contains a single BCL2 homology domain 3 (BH3), and has been shown to bind BCL2 proteins and function as an apoptotic activator. This protein is found to be sequestered to myosin V motors by its association with dynein light chain 2, which may be important for sensing intracellular damage and triggering apoptosis. Alternatively spliced transcript variants encoding different isoforms have been identified.

Recommended Dilutions

Immunogen Information

Gene ID
90427

Swiss Prot
Q96LC9

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-184 of human BMF (NP_277038.1).

Synonyms
BMF

Contact

 | www.abclonal.com

Product Information

Source
Rabbit

Isotype
IgG

Purification
Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.