

A1310

Leader in Biomolecular Solutions for Life Science



HFE Rabbit pAb

Catalog No.: A1310

1 Publications

Basic Information

Observed MW

50kDa

Calculated MW

40kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene is a membrane protein that is similar to MHC class I-type proteins and associates with beta2-microglobulin (beta2M). It is thought that this protein functions to regulate iron absorption by regulating the interaction of the transferrin receptor with transferrin. The iron storage disorder, hereditary haemochromatosis, is a recessive genetic disorder that results from defects in this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3077

Swiss Prot

Q30201

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 80-306 of human HFE (NP_000401.1).

Synonyms

HH; HFE1; HLA-H; MVCD7; TFQTL2; HFE

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using HFE antibody (A1310) at 1:1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.