

MORN1 Rabbit pAb

Catalog No.: A13194

Basic Information

Observed MW

56kDa

Calculated MW

54kDa

Category

Primary antibody

Applications

ELISA, WB

Cross-Reactivity

Human

Background

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79906

Swiss Prot

Q5T089

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 201-350 of human MORN1 (NP_001287989.1).

Synonyms

MORN1

Contact

 | www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of extracts of K-562 cells, using MORN1 antibody (A13194) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.