

A13263

Leader in Biomolecular Solutions for Life Science



VGLL3 Rabbit pAb

Catalog No.: A13263

Basic Information

Observed MW

36kDa

Calculated MW

36kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Rat

Background

Predicted to enable protein C-terminus binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to be active in nucleus.

Recommended Dilutions

WB 1:1000 - 1:5000

Immunogen Information

Gene ID

389136

Swiss Prot

A8MV65

Immunogen

A synthetic peptide corresponding to a sequence within amino acids 226-326 of human VGLL3 (NP_057290.2).

Synonyms

VGL3; VGL-3; VGLL3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates, using VGLL3 Rabbit pAb (A13263) at 1:2000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25ug per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 30s.