

A13413

Leader in Biomolecular Solutions for Life Science



LHX6 Rabbit pAb

Catalog No.: A13413

Basic Information

Observed MW

Calculated MW

40kDa

Category

Mouse Monoclonal Antibody

Applications

ELISA

Cross-Reactivity

Human

Background

This gene encodes a member of a large protein family that contains the LIM domain, a unique cysteine-rich zinc-binding domain. The encoded protein has tandem LIM domains as well as a DNA-binding homeodomain. The protein functions as a transcription factor involved in embryogenesis and head development and is highly expressed in neural crest derived mesenchyme cells. Alternative splicing results in multiple transcript variants encoding distinct isoforms.

Recommended Dilutions

Immunogen Information

Gene ID

26468

Swiss Prot

Q9UPM6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 143-392 of human LHX6 (NP_055183.2).

Synonyms

hLHX6; LHX6.1; LHX6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.