

A13545

Leader in Biomolecular Solutions for Life Science



UAP1 Rabbit pAb

Catalog No.: A13545

Basic Information

Observed MW

Calculated MW

59kDa

Category

Polyclonal Antibody

Applications

IHC-P,ELISA

Cross-Reactivity

Human,Rat

Background

Enables identical protein binding activity. Predicted to be involved in UDP-N-acetylglucosamine biosynthetic process. Located in cytosol; nucleoplasm; and plasma membrane.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

6675

Swiss Prot

Q16222

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-290 of human UAP1 (NP_003106.3).

Synonyms

AGX; AGX1; AGX2; SPAG2; UAP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunohistochemistry analysis of UAP1 in paraffin-embedded rat brain using UAP1 Rabbit pAb (A13545) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.



Immunohistochemistry analysis of UAP1 in paraffin-embedded human gastric cancer using UAP1 Rabbit pAb (A13545) at dilution of 1:200 (40x lens). Perform microwave antigen retrieval with 10 mM PBS buffer pH 7.2 before commencing with IHC staining protocol.