

A13604

Leader in Biomolecular Solutions for Life Science



RFFL Rabbit pAb

Catalog No.: A13604

Basic Information

Observed MW

Calculated MW
41kDa

Category
Polyclonal Antibody

Applications
ELISA

Cross-Reactivity
Human

Background

Enables enzyme binding activity; p53 binding activity; and ubiquitin protein ligase activity. Involved in cellular protein metabolic process; negative regulation of cysteine-type endopeptidase activity involved in execution phase of apoptosis; and negative regulation of signal transduction. Located in endosome membrane and plasma membrane.

Recommended Dilutions

Immunogen Information

Gene ID	Swiss Prot
117584	Q8WZ73

Immunogen
Recombinant fusion protein containing a sequence corresponding to amino acids 1-190 of human RFFL (NP_001017368.1).

Synonyms
CARP2; FRING; CARP-2; RNF189; RNF34L; RIFIFYLIN; RFFL

Contact

 | www.abclonal.com

Product Information

Source	Isotype	Purification
Rabbit	IgG	Affinity purification

Storage
Store at -20°C. Avoid freeze / thaw cycles.
Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.