

A13760

Leader in Biomolecular Solutions for Life Science



IFIT2 Rabbit pAb

Catalog No.: A13760

1 Publications

Basic Information

Observed MW

50kDa

Calculated MW

55kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Enables RNA binding activity. Involved in negative regulation of protein binding activity; positive regulation of apoptotic process; and response to virus. Located in endoplasmic reticulum.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

3433

Swiss Prot

P09913

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 273-472 of human IFIT2 (NP_001538.4).

Synonyms

P54; G10P2; IFI54; ISG54; cig42; IFI-54; IFIT-2; ISG-54; GARG-39; IFI-54K; ISG-54K; ISG-54 K; IFIT2

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of various lysates using IFIT2 Rabbit pAb (A13760) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.