

A13774

Leader in Biomolecular Solutions for Life Science



JAKMIP1 Rabbit pAb

Catalog No.: A13774

Basic Information

Observed MW

73kDa

Calculated MW

73kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

Enables GABA receptor binding activity and RNA binding activity. Involved in cognition. Is extrinsic component of membrane. Part of ribonucleoprotein complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

152789

Swiss Prot

Q96N16

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 337-626 of human JAKMIP1 (NP_653321.1).

Synonyms

JAMIP1; MARLIN1; Gababrbp; JAKMIP1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using JAKMIP1 Rabbit pAb (A13774) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (A5014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 90s.