

# ZNF227 Rabbit pAb

Catalog No.: A13793

## Basic Information

### Observed MW

Refer to figures

### Calculated MW

92kDa

### Category

Primary antibody

### Applications

ELISA, WB

### Cross-Reactivity

Human

## Background

Predicted to enable DNA-binding transcription activator activity, RNA polymerase II-specific and RNA polymerase II cis-regulatory region sequence-specific DNA binding activity. Predicted to be involved in regulation of transcription by RNA polymerase II. Predicted to act upstream of or within negative regulation of transcription, DNA-templated. Predicted to be active in nucleus.

## Recommended Dilutions

**WB** 1:500 - 1:2000

## Immunogen Information

### Gene ID

7770

### Swiss Prot

Q86WZ6

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 60-230 of human ZNF227 (NP\_872296.1).

### Synonyms

ZNF227

## Contact

 | [www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification

### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.