

A13812

Leader in Biomolecular Solutions for Life Science



ZMYM4 Rabbit pAb

Catalog No.: A13812

Basic Information

Observed MW

240kDa

Calculated MW

173kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable DNA binding activity. Involved in cytoskeleton organization and regulation of cell morphogenesis.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9202

Swiss Prot

Q5VZL5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 35-230 of human ZMYM4 (NP_005086.2).

Synonyms

MYM; CDIR; ZNF262; ZNF198L3; ZMYM4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using ZMYM4 Rabbit pAb (A13812) at 1:3000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 90s.