

A13982

Leader in Biomolecular Solutions for Life Science



FSHR Rabbit pAb

Catalog No.: A13982

Basic Information

Observed MW

Refer to figures

Calculated MW

78kDa

Category

Mouse Monoclonal Antibody

Applications

WB

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene belongs to family 1 of G-protein coupled receptors. It is the receptor for follicle stimulating hormone and functions in gonad development. Mutations in this gene cause ovarian dysgenesis type 1, and also ovarian hyperstimulation syndrome. Alternative splicing results in multiple transcript variants.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

2492

Swiss Prot

P23945

Immunogen

Recombinant protein of human FSHR

Synonyms

LGR1; ODG1; FSHR1; FSHRO; FSHR

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.