

A14034

Leader in Biomolecular Solutions for Life Science



OMG Rabbit pAb

Catalog No.: A14034

Basic Information

Observed MW

Refer to figures

Calculated MW

50kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to enable identical protein binding activity. Predicted to be involved in neuron projection regeneration. Predicted to act upstream of or within regulation of collateral sprouting of intact axon in response to injury. Predicted to be located in plasma membrane.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4974

Swiss Prot

P23515

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 190-420 of human OMG (NP_002535.3).

Synonyms

OMGP; OMG

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.