

A14088

Leader in Biomolecular Solutions for Life Science



NUMBL Rabbit pAb

Catalog No.: A14088

Basic Information

Observed MW

Calculated MW

65kDa

Category

Mouse Monoclonal Antibody

Applications

IHC-P,ELISA

Cross-Reactivity

Human

Background

Involved in cytokine-mediated signaling pathway and protein metabolic process. Predicted to be active in cytoplasm.

Recommended Dilutions

IHC-P 1:50 - 1:200

Immunogen Information

Gene ID

9253

Swiss Prot

Q9Y6R0

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 330-609 of human NUMBL (NP_004747.1).

Synonyms

NBL; CAG3A; CTG3a; NUMBR; NUMB-R; TNRC23; NUMBLIKE; NUMBL

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.