

A14148

Leader in Biomolecular Solutions for Life Science



BCL2L14/BCLG Rabbit pAb

Catalog No.: A14148

Basic Information

Observed MW

Refer to figures

Calculated MW

37kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

The protein encoded by this gene belongs to the BCL2 protein family. BCL2 family members form hetero- or homodimers and act as anti- or pro-apoptotic regulators that are involved in a wide variety of cellular activities. Overexpression of this gene has been shown to induce apoptosis in cells. Three alternatively spliced transcript variants encoding two distinct isoforms have been reported for this gene.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

79370

Swiss Prot

Q9BZR8

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-327 of human BCL2L14/BCLG (NP_620049.1).

Synonyms

BCLG; BCL2L14/BCLG

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.