

A14196

Leader in Biomolecular Solutions for Life Science



## WNK1 Rabbit pAb

Catalog No.: A14196

### Basic Information

**Observed MW**

300kDa

**Calculated MW**

251kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Mouse, Rat

### Background

This gene encodes a member of the WNK subfamily of serine/threonine protein kinases. The encoded protein may be a key regulator of blood pressure by controlling the transport of sodium and chloride ions. Mutations in this gene have been associated with pseudohypoaldosteronism type II and hereditary sensory neuropathy type II. Alternatively spliced transcript variants encoding different isoforms have been described but the full-length nature of all of them has yet to be determined.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

65125

**Swiss Prot**

Q9H4A3

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 1-210 of human WNK1 (NP\_055638.2).

**Synonyms**

KDP; PSK; p65; HSN2; HSN2; PRKWKN1; PPP1R167; WNK1

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.

Validation Data



Western blot analysis of extracts of various cell lines, using WNK1 antibody (A14196) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 90s.