

A14418

Leader in Biomolecular Solutions for Life Science



## NLGN2 Rabbit pAb

Catalog No.: A14418

### Basic Information

#### Observed MW

Refer to Figures

#### Calculated MW

91kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Mouse, Rat

### Background

This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

57555

#### Swiss Prot

Q8NFZ4

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 15-180 of human NLGN2 (NP\_065846.1).

#### Synonyms

NLGN2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.