

A14488

Leader in Biomolecular Solutions for Life Science



NOV Rabbit pAb

Catalog No.: A14488

Basic Information

Observed MW

Refer to Figures

Calculated MW

39kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Mouse

Background

The protein encoded by this gene is a small secreted cysteine-rich protein and a member of the CCN family of regulatory proteins. CCN family proteins associate with the extracellular matrix and play an important role in cardiovascular and skeletal development, fibrosis and cancer development.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4856

Swiss Prot

P48745

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 71-357 of human NOV (NP_002505.1).

Synonyms

NOV; NOVh; IBP-9; IGFBP9; IGFBP-9

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.