

A14518

Leader in Biomolecular Solutions for Life Science



## POMGNT2 Rabbit pAb

Catalog No.: A14518

### Basic Information

**Observed MW**

75kDa

**Calculated MW**

67kDa

**Category**

Polyclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Human, Mouse, Rat

### Background

This gene encodes a protein with glycosyltransferase activity although its function is not currently known.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

84892

**Swiss Prot**

Q8NAT1

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 160-340 of human POMGNT2 (NP\_116195.2).

**Synonyms**

AGO61; GTDC2; MDDGA8; MDDGC8; C3orf39; POMGNT2

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**


Affinity purification

**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

# Validation Data



Western blot analysis of various lysates using POMGNT2 Rabbit pAb (A14518) at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution. Lysates/proteins: 25µg per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit (RM00020). Exposure time: 1s.