

A15029

Leader in Biomolecular Solutions for Life Science



## SLC25A6 Rabbit pAb

Catalog No.: A15029

### Basic Information

#### Observed MW

#### Calculated MW

33kDa

#### Category

Polyclonal Antibody

#### Applications

ELISA

#### Cross-Reactivity

Human

### Background

This gene is a member of the mitochondrial carrier subfamily of solute carrier protein genes. The product of this gene functions as a gated pore that translocates ADP from the cytoplasm into the mitochondrial matrix and ATP from the mitochondrial matrix into the cytoplasm. The protein is implicated in the function of the permeability transition pore complex (PTPC), which regulates the release of mitochondrial products that induce apoptosis. The human genome contains several non-transcribed pseudogenes of this gene.

### Recommended Dilutions

### Immunogen Information

#### Gene ID

293

#### Swiss Prot

P12236

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-298 of human SLC25A6 (NP\_001627.2).

#### Synonyms

ANT; AAC3; ANT3; ANT 2; ANT 3; ANT3Y; SLC25A6

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.