

A15038

Leader in Biomolecular Solutions for Life Science



CHD1 Rabbit pAb

Catalog No.: A15038

Basic Information

Observed MW

Refer to figures

Calculated MW

197kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

The CHD family of proteins is characterized by the presence of chromo (chromatin organization modifier) domains and SNF2-related helicase/ATPase domains. CHD genes alter gene expression possibly by modification of chromatin structure thus altering access of the transcriptional apparatus to its chromosomal DNA template.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

1105

Swiss Prot

O14646

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1501-1710 of human CHD1 (NP_001261.2).

Synonyms

CHD-1; PILBOS; CHD1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.