

A15081

Leader in Biomolecular Solutions for Life Science



PLD1 Rabbit pAb

Catalog No.: A15081

1 Publications

Basic Information

Observed MW

Refer to figures

Calculated MW

124kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

This gene encodes a phosphatidylcholine-specific phospholipase which catalyzes the hydrolysis of phosphatidylcholine in order to yield phosphatidic acid and choline. The enzyme may play a role in signal transduction and subcellular trafficking. Alternative splicing results in multiple transcript variants with both catalytic and regulatory properties.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5337

Swiss Prot

Q13393

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 875-1036 of human PLD1 (NP_001123553.1).

Synonyms

CVDD; CVDP1; PLD1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.