

A15503

Leader in Biomolecular Solutions for Life Science



XPO4 Rabbit pAb

Catalog No.: A15503

Basic Information

Observed MW

Refer to figures

Calculated MW

130kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

XPO4 belongs to a large family of karyopherins (see MIM 602738) that mediate the transport of proteins and other cargo between the nuclear and cytoplasmic compartments (Lipowsky et al., 2000 [PubMed 10944119]).

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

64328

Swiss Prot

Q9C0E2

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 400-700 of human XPO4 (NP_071904.4).

Synonyms

exp4; XPO4

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.