

A15521

Leader in Biomolecular Solutions for Life Science



UBXN6 Rabbit pAb

Catalog No.: A15521

Basic Information

Observed MW

Refer to figures

Calculated MW

50kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Involved in ERAD pathway; endosome to lysosome transport via multivesicular body sorting pathway; and macroautophagy. Located in bounding membrane of organelle and cytosol. Is extrinsic component of membrane. Part of endosome and protein-containing complex.

Recommended Dilutions

WB 1:200 - 1:2000

Immunogen Information

Gene ID

80700

Swiss Prot

Q9BZV1

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 49-388 of human UBXN6 (NP_001164562.1).

Synonyms

UBXD1; UBXDC2; UBXN6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.