

A15648

Leader in Biomolecular Solutions for Life Science



BMP6 Rabbit pAb

Catalog No.: A15648

Basic Information

Observed MW

Refer to figures

Calculated MW

57kDa

Category

Mouse Monoclonal Antibody

Applications

IF/ICC,ELISA

Cross-Reactivity

Human,Mouse

Background

This gene encodes a secreted ligand of the TGF-beta (transforming growth factor-beta) superfamily of proteins. Ligands of this family bind various TGF-beta receptors leading to recruitment and activation of SMAD family transcription factors that regulate gene expression. The encoded preproprotein is proteolytically processed to generate each subunit of the disulfide-linked homodimer. This protein regulates a wide range of biological processes including iron homeostasis, fat and bone development, and ovulation. Differential expression of this gene may be associated with progression of breast and prostate cancer. Mutations in this gene may be associated with iron overload in human patients.

Recommended Dilutions

IF/ICC 1:100 - 1:500

Immunogen Information

Gene ID

654

Swiss Prot

P22004

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 375-513 of human BMP6 (NP_001709.1).

Synonyms

IO; VGR; VGR1; BMP6

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

Validation Data



Immunofluorescence analysis of U2OS cells using BMP6 Rabbit pAb (A15648) at dilution of 1:100. Secondary antibody: Cy3 Goat Anti-Rabbit IgG (H+L) (AS007) at 1:500 dilution. Blue: DAPI for nuclear staining.