

A15702

Leader in Biomolecular Solutions for Life Science



Kallistatin (SERPINA4) Rabbit pAb

Catalog No.: A15702

Basic Information

Observed MW

Refer to figures

Calculated MW

49kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Predicted to enable serine-type endopeptidase inhibitor activity. Predicted to be involved in negative regulation of endopeptidase activity. Located in extracellular exosome. Biomarker of diabetic retinopathy.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

5267

Swiss Prot

P29622

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 178-427 of human Kallistatin (SERPINA4) (NP_006206.2).

Synonyms

KAL; KST; PI4; KLST; PI-4; kallistatin; Kallistatin (SERPINA4)

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide, 50% glycerol, pH 7.3.