

A15806

Leader in Biomolecular Solutions for Life Science



## TAZ Rabbit pAb

Catalog No.: A15806

2 Publications

### Basic Information

#### Observed MW

55kDa

#### Calculated MW

44kDa

#### Category

Mouse Monoclonal Antibody

#### Applications

WB,IF/ICC,ELISA

#### Cross-Reactivity

Human,Mouse

### Background

Enables transcription coactivator activity. Involved in several processes, including hippo signaling; positive regulation of cell differentiation; and regulation of signal transduction. Located in cytosol and nuclear body.

### Recommended Dilutions

WB	1:500 - 1:2000
IF/ICC	1:50 - 1:200

### Immunogen Information

#### Gene ID

25937

#### Swiss Prot

Q9GZV5

#### Immunogen

A synthetic peptide corresponding to a sequence within amino acids 300-400 of human TAZ (NP\_056287.1).

#### Synonyms

TAZ

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification


Affinity purification

#### Storage


Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.

## Validation Data



Western blot analysis of extracts of various cell lines, using TAZ antibody (A15806) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 30s.



Immunofluorescence analysis of U2OS cells using TAZ antibody (A15806) at dilution of 1:100. Blue: DAPI for nuclear staining.