

A15927

Leader in Biomolecular Solutions for Life Science



C12orf29 Rabbit pAb

Catalog No.: A15927

Basic Information

Observed MW

37kDa

Calculated MW

37kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse

Background

Predicted to act upstream of or within hematopoietic progenitor cell differentiation.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

91298

Swiss Prot

Q8N999

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-325 of human C12orf29 (NP_001009894.2).

Synonyms

C12orf29

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.

Validation Data



Western blot analysis of various lysates using C12orf29 Rabbit pAb (A15927) at 1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25µg per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Basic Kit (RM00020).
Exposure time: 10s.