

A15933

Leader in Biomolecular Solutions for Life Science



MVB12A Rabbit pAb

Catalog No.: A15933

Basic Information

Observed MW

29kDa

Calculated MW

29kDa

Category

Polyclonal Antibody

Applications

WB,ELISA

Cross-Reactivity

Rat

Background

Enables lipid binding activity and ubiquitin binding activity. Involved in regulation of epidermal growth factor receptor signaling pathway; viral budding; and virus maturation. Located in several cellular components, including Golgi apparatus; centrosome; and nucleoplasm. Part of ESCRT I complex.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

93343

Swiss Prot

Q96EY5

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-273 of human MVB12A (NP_612410.1).

Synonyms

CFBP; FAM125A; MVB12A

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification


Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

Validation Data



Western blot analysis of lysates from rat pancreas, using MVB12A Rabbit pAb (A15933) at 1000 dilution.
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.
Lysates/proteins: 25 μ g per lane.
Blocking buffer: 3% nonfat dry milk in TBST.
Detection: ECL Enhanced Kit (RM00021).
Exposure time: 10s.