

A16000

Leader in Biomolecular Solutions for Life Science



# NUSAP1 Rabbit pAb

Catalog No.: A16000

1 Publications

## Basic Information

### Observed MW

49kDa

### Calculated MW

49kDa

### Category

Polyclonal Antibody

### Applications

WB,ELISA

### Cross-Reactivity

Human,Mouse

## Background

NUSAP1 is a nucleolar-spindle-associated protein that plays a role in spindle microtubule organization (Raemaekers et al., 2003 [PubMed 12963707]).

## Recommended Dilutions

WB 1:500 - 1:2000

## Immunogen Information

### Gene ID

51203

### Swiss Prot

Q9BXS6

### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-61 of human NUSAP1 (NP\_060924.4).

### Synonyms

LNP; ANKT; SAPL; BM037; NUSAP; Q0310; PRO0310p1; NUSAP1

## Contact



[www.abclonal.com](http://www.abclonal.com)

## Product Information

### Source

Rabbit

### Isotype

IgG

### Purification

Affinity purification


### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal,50% glycerol,pH7.3.

## Validation Data

---



Western blot analysis of various lysates using NUSAP1 Rabbit pAb (A16000) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Basic Kit (RM00020).  
Exposure time: 60s.