

A16031

Leader in Biomolecular Solutions for Life Science



## C1QTNF7 Rabbit pAb

Catalog No.: A16031

### Basic Information

#### Observed MW

Refer to figures

#### Calculated MW

31kDa

#### Category

Polyclonal Antibody

#### Applications

WB, ELISA

#### Cross-Reactivity

Mouse, Rat

### Background

Predicted to enable identical protein binding activity. Predicted to be located in extracellular space. Predicted to be part of collagen trimer.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

#### Gene ID

114905

#### Swiss Prot

Q9BXJ2

#### Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 24-190 of human C1QTNF7 (NP\_001128642.1).

#### Synonyms

CTRP7; ZACRP7; C1QTNF7

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

#### Source

Rabbit

#### Isotype

IgG

#### Purification

Affinity purification

#### Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH7.3.