

A16170

Leader in Biomolecular Solutions for Life Science



## EFCAB3 Rabbit pAb

Catalog No.: A16170

### Basic Information

**Observed MW**

55kDa

**Calculated MW**

50kDa

**Category**

Mouse Monoclonal Antibody

**Applications**

WB, ELISA

**Cross-Reactivity**

Mouse, Rat

### Background

Predicted to enable calcium ion binding activity.

### Recommended Dilutions

WB 1:500 - 1:2000

### Immunogen Information

**Gene ID**

146779

**Swiss Prot**

Q8N7B9

**Immunogen**

Recombinant fusion protein containing a sequence corresponding to amino acids 249-438 of human EFCAB3 (NP\_775774.1).

**Synonyms**

EFCAB3

### Contact



[www.abclonal.com](http://www.abclonal.com)

### Product Information

**Source**

Rabbit

**Isotype**

IgG

**Purification**

Affinity purification


**Storage**

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.

## Validation Data

---



Western blot analysis of various lysates using EFCAB3 Rabbit pAb (A16170) at 1:1000 dilution.  
Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) (AS014) at 1:10000 dilution.  
Lysates/proteins: 25µg per lane.  
Blocking buffer: 3% nonfat dry milk in TBST.  
Detection: ECL Enhanced Kit (RM00021).  
Exposure time: 90s.