

A16186

Leader in Biomolecular Solutions for Life Science



XPR1 Rabbit pAb

Catalog No.: A16186

Basic Information

Observed MW

Refer to figures

Calculated MW

82kDa

Category

Polyclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human, Mouse, Rat

Background

The protein encoded by this gene is a receptor for the xenotropic and polytropic classes of murine leukemia viruses. The encoded protein is involved in phosphate homeostasis by mediating phosphate export from the cell. Defects in this gene have been associated with idiopathic basal ganglia calcification-6.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

9213

Swiss Prot

Q9UBH6

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 464-631 of human XPR1 (NP_001129141.1).

Synonyms

X3; SYG1; IBGC6; SLC53A1; XPR1

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.