

A16238

Leader in Biomolecular Solutions for Life Science



NDUFA3 Rabbit pAb

Catalog No.: A16238

Basic Information

Observed MW

Refer to figures

Calculated MW

9kDa

Category

Mouse Monoclonal Antibody

Applications

WB, ELISA

Cross-Reactivity

Human

Background

Involved in mitochondrial respiratory chain complex I assembly. Located in mitochondrion. Part of mitochondrial respiratory chain complex I.

Recommended Dilutions

WB 1:500 - 1:2000

Immunogen Information

Gene ID

4696

Swiss Prot

O95167

Immunogen

Recombinant fusion protein containing a sequence corresponding to amino acids 1-84 of human NDUFA3 (NP_004533.1).

Synonyms

B9; CI-B9; NDUFA3

Contact



www.abclonal.com

Product Information

Source

Rabbit

Isotype

IgG

Purification

Affinity purification

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Buffer: PBS with 0.01% thimerosal, 50% glycerol, pH 7.3.